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The straight and 30° angled IDUnit abutment is available for all the IDI range with different gingival heights.

The IDUnit components are dedicated to the IDI implants range only.

The I.D.I (Implants Diffusion International) firm accepts no liability for the use of products manufactured by a third part.

It is the clinician's responsibility to determinate the correct use of the IDUnit and its application.

Implants Diffusion International can not be responsible for direct and/or indirect damages following a wrong appreciation with the use of the IDUnit connections.

The IDI firm (Implants Diffusion International) accepts no liability for the damages caused by the use of the IDUnit.



The IDUnit abutments are intended for the fixed prostheses and for the screw-on removable prostheses.

It is recommended to attend a training on the use of the IDUnit abutments.

The success of a prosthetic realization depends on the collaboration between the practitioners and the dental technician.

It is recommended to master the occlusion and the distribution of the occlusal stress for the realization of the prostheses on the IDUnit.

The screw-on metallic frameworks on the IDUnit abutments must be rigid enough to oppose to any flexions and to any deformations linked to the ceramics heating up.

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Elements constituting the IDUnit concept:

- A) Straight IDUnit abutment;
- B) Angled 30° IDUnit abutment;
- C) Rotational impression cylinder for IDUnit abutment (Ref.321);
- D) Protection cap for IDUnit abutment (Ref.335);
- E) IDUnit analogue (Ref.333);
- F) IDUnit burnout cylinder (Ref.336);
- G) IDUnit retaining screw (Ref.216).

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Method of prosthetic realization:

Indirect method

- 1-Take the impression in mouth with the rotationally secure impression cylinder on the implant, close or open-tray technique (Ref. H904, H904L, 2004, 2004F).
- 2-The laboratory positions the implant analogue onto the impression cylinder, it casts the artificial gum and prepares the model in plaster.
- 3- The technician screws the IDUnit abutments on the implant analogues.
- 4-The technician prepares the screw-on bar thanks to the burnout cylinders (F).
- 5-The technician prepares the cosmetic element on the bar, in acrylics, in ceramic-composite, or in ceramics.
- 6- In the case of an angled IDUnit is used, the technician will specify the exact position of the angled IDUnit on the implant in the model, and if need, he will realize a repositioning key.

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Direct method

- 1-The practitioner screws in mouth the IDUnit abutments (A or B) on the implants to 25 N.cm, with the screwdriver reference 0114 or 0014 or 1014 or 1114 so that the retaining screws are positioned lingually as regards the prosthetic elements.
 - 2-The practitioner screws the impression cylinder (C) to 5 N.cm with the same screwdriver as previously mentioned. He fills the upper part of the impression cylinder with wax and takes the impression with high quality elastomeric material.
 - 3-The practitioner unscrews the impression cylinders and repositions them onto the silicone impression. Then, he sends them to the technician.
 - 4-The practitioner fixes the protection caps (D) on the IDUnit abutments fixed in mouth.
 - 5-The technician screws the IDUnit abutment analogues (E) onto the impression cylinders, he casts the artificial silicone gum and prepares the model in plaster.
 - 6-The technician realizes the screw-on bar thanks to the burnout cylinders (F).
 - 7-The technician realizes the cosmetic element on the bar in acrylics, ceramo-composite or in ceramics.
- The practitioner screws in mouth the screw-on prosthesis (G) with the short screws to 15 N.cm on the IDUnit abutments after having removed the protection caps.